LSM 880 Basic Operation



MultiChannel & TimeSeries & Z-Stack & TileScan



2016-05-15



1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	



1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	



Startup of the System





1、打开稳压电源:

- a、打开稳压电源背后总开关;
- b、打开稳压电源前面开关;
- c、确定稳压电源显示220V;



Startup of the System



- 1 Main switch ON/OFF
- 2 COMPONENTS switch ON/OFF
- **3** SYSTEM PC switch ON/OFF
- 4 Key switch

*在打开整个系统电源之前,确保 钥匙(图中4)处于"ON";
2、打开"MAIN SWITCH"(图1);
3、打开"SYSTEM/PC"(图3),
随后打开电脑,进入系统;
4、打开"COMPONENTS"(图2);
5、打开金属卤化物灯;



ZEISS

Startup of the System



7、点击"Start System"



Startup of the System





- 8、如果需要使用458、488、514
 需要打开Argon激光器:

 a、在"Laser"中"Argon"栏
 选择"Standby",等待激光器预热;
 b、预热结束后"status"显示
 "Ready"后选择"On";
- 9、随后可以开始共聚焦成像



Shutdown of the System

A reset		۵ ا
Laser	Laser Lines [nm]	Power
À Argon	458, 488, 514	On 🔻
À Diode 405-30	405	
📤 DPSS 561-10	561	On 🔻
HeNe633	633	Off
Laser Properties	458, 488, 514 nm	
Laser Properties Wavelength Status	458, 488, 514 nm Ready	
Laser Properties Wavelength Status Tube Current	458, 488, 514 nm Ready 5.2 A	

关机过程基本类似开机的倒序过程:

- 1、关闭金属卤化物灯;
- 2、关闭激光器,在 "power" 中选择 "off",关闭已打开的激光器:
- 3、关闭ZEN软件;
- 4、关闭电脑;
- 5、关闭电脑后等待激光器冷却约5min, 此时Argon激光器的风扇停转;
- 6、按顺序关闭"COMPONENTS"→
 - "SYSTEM/PC" → "MAIN SWITCH"
- 7、关闭稳压电源。



1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	

Acquiring Multi-Channel images



ZEINS

Acquiring Multi-Channel images

(IIN) ZEN 2.1	
File View Acquisition Maintain Macro Tools Window Help	
Locate Acministran	Image 1 🛛
Experiment Manager	2D
Smart Setup Vision Visi	2.5D
AF Image: Continuous Image: Continuous <t< th=""><th>Histo</th></t<>	Histo
I Z-Stack □ Time Series □ Bleaching □ Tile Scan □ Tile Scan 256 kB 256 kB Speed	Prohie Linfo
Positions Averaging Regions → Start Experiment Averaging Number 2 → Bit Depth 8 Bit →	

2、进入 "Acquisition" 界面; 3、选择之前保存好的光路设置→ "channel"和 "Acquisition mode"设置 4、或者新建光路设置→ "Smart Setup"

Switch track every Line	Pixel Size: 0.42 µm ↔ 1 0.0 ℃ C ↓ 0.0 ℃ C	
	Zoom 1 10 1 Reset All	
400 500 600 700		
Use Dye Color Detector Range +	✓ A Channels ✓ Show all	AAA
□ Alexa Fluor 488 ▼ □ ▼ Ch1 493-630nm▼ -		
□ ▼ ■ Ch2 GaAsP 415-735nm▼	Tracks Channels	
□ ▼ □ ▼ Ch3 415-735nm▼	Track 1 T PMT	
Reflection		
MBS 488 🕖 Visible Light 📼		Show all
Plate 🕖 Invisible Light 🍊	Track Configuration not defined	Zoom 101% Q Q 153 ÷
None	Track 1 - LSM	Tools Q□ Q; → VInterpolation
ZEISS, Lu XI, MIC		2016-05-15

ZEINN



Channels





6、在"live"下设置Channels中的激光强度"Laser",针孔大小"Pinhole",检测器"Gain"值,以及"digital gain"或"digital offset";每个track单独设置,所以每次只勾选一个track,并选中该track(选中track高亮);
A、"Pinhole"一般设置为1 AU,增大Pinhole可以提高图像亮度,但会增加非焦面信息;减少Pinhole可以增加景深,但是会减少图像亮度;
B、"Gain"和"Digital Gain"增加可以

B、"Gain"和"Digital Gain"增加可以 增加图像亮度,但是也会提高背景噪音。









Acquisition Mode



🔹 🛥 Acquisition Mode		Show all 📝
Objective Plan-	Apochromat 10x/0.45 M27	
Scan Mode Fram	e 🔻	
Frame Size X 5	12 🗘 X * Y	Y 512
Line Step 1		Optimal
Speed ——	7	‡) Max
Pixel Dwell 2.06 µse	c Scan Time 1.27 sec	
Averaging		
Number 2	Bit Depth	8 Bit
Mode Line	Direction	<> v
Method Mea	n 🔽 Corr X ——	0.00
	Corr Y —— 🗍	0.00
		Auto
HDR		
Airyscan Multitrack		ILEX Setup
😇 scan Area		
	Image Size: 848.5 μm	x 848.5 μm
• • •	Pixel Size: 1.66 μm	
		0.0 🗘 C
		0.0 🗘 C
	φ — <u></u>	0.0 🗘 🛈
	Zoom () [1	.0 🗘 1
		Reset All
	Line Select	

7、在Acquisition Mode下主要设置如下参数: A、通过scan area选择扫描区域或通过图像窗 口下的"crop"选择扫描区域;

Dimensions	Display Graphics		
Zoom		━,	76 🗘
Tools	Q Q_‡ 📑		✓ Interpolation
Channels	Ch1 T PMT		
	Single Channel	✓ Range Indicator	Quick Color Setup
✓ Merged		Reuse 🎒 Croj	p Positions Stage

Acquisition Mode



👻 🛥 Acquisition M	lode	Show all
Objective	Plan-Apochromat 10x/0.45 M2	7
Scan Mode	Frame	
	X 512 🗘 X*Y	Y 512
Line Step	1	Optimal
	7	‡ Max
Pixel Dwell 2.	.06 µsec Scan Time 1.27 sec	
Averaging		
Number	2 Bit Depth	8 Bit 🔻
Mode	Line Direction	<> •
Method	Mean Corr X -	0.00 🗘
	Corr Y —— 🗍	0.00
		Auto
HDR		
Airyscan Multitrack		ILEX Setup
Scan Area		
<u></u>	Image Size: 848.5 μm	x 848.5 μm
	Pixel Size: 1.66 μm	
	200m []	Reset All
C	Line Select	

7F

- B、设置speed:扫描速度越慢,信噪比越好, 但光漂白越多;
- C、Averaging: 增加averaging次数可以减少 噪音,但会增加扫描时间;
- D、Direction: 双向扫描可以减少扫描时间;
- E、Frame Size: 一般选择512×512或
 - 1024×1024,图像越大,扫描时间越长

Acquiring Multi-Channel images





Experiment Manager and Reuse





可以通过Experiment Manager来保存拍摄的参数,或打开已经保存的图片(czi格式),通过"Reuse"来调用上次拍摄的参数设置



1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	

First/Last



- 1、选择Z-stack;
- 2、在live下调节焦距选择层扫图像的上下范围: "set first" 和 "set last";
- 3、单击"optimal"让间距在最适合范围;
- 4、单击"Start Experiment"



Center



 1、center模式下, live下选择成 像的中间位置,单击"center"
 ,然后设置需要层扫的层数
 Slices,并单击"optimal";
 2、单击"Start Experiment"



Match Pinhole



多通道荧光拍摄Z-stack需 要考虑光切厚度不一致的 问题:

1、可以通过点击"Match Pinhole"自动调节不同 track的针孔使光切厚度相 似;

*这种方法的缺点在于可能 会使长波长的针孔过于小 ,不利于弱荧光成像。





Match Pinhole



2、通过手动调节针孔 到一致,可以保证荧 光强度的同时,保证 光切厚度一致。







1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	

Time Series image





- 1、选择"Time Series";
- 2、拍摄要进行的循环数cycles;
- 3、循环之间的间隔Interval;
 - interval:两次循环开始时间的间隔,因此interval包含了上一个循环的拍摄时间;
- 4, "Start Experiment"



1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	

Tile Scan





- 1、选择"Tile Scan";
- 2、Centered grid下以现在视野为中心,拍摄拼图:
- A、Horizontal & Vertical 水平和垂直拍摄范围;
- B、拼图之间的重叠overlap、是否双向拼图Bi-directional以及拍摄后拼接 Online stitching;
- 3、 "Start Experiment"

Tile Scan



Bounding grid

👻 🎟 Tile Scan 🗸 Show all 📝				
Center	ed grid	Bounding grid	Con	vex hull
	Tiles	Pixels	Size	
Horizontal	11	5120	749.73 μm	
Vertical		1894	277.34 μm	
Overlap	10 ‡	% Add	Remove	Remove all
Rotation	0.0000 ‡			
Bi-directional Online stitching				
Number	x [μm]	y [µm]		
1	-7542.3	-1806	.440	
2	-7752.9 -82 <u>06.</u> 7	764 -1970 764 -1 <u>970</u>	.115	
			Load	Save
Scan overview image				

Bounding grid下通过移动载物台到想拼 接的图像边缘,单击"Add",添加边缘范围。 软件自动定义拼图范围。 "Start Experiment"。

Tile Scan



Convex hull

👻 🎫 Tile Scan 🛛 🖍 Show all 📝				
Centere	ed grid	Bounding grid		Convex hull
	Tiles	Pixels	Size	
Horizontal	11	5120	749.73	μm
Vertical		1894	277.34	μm
Overlap	10 ‡	% Add	Remo	ove Remove all
Rotation	0.0000 ‡			
✓ Bi-directi Online st	ional titching			
😌 Marked	positions			
Number	x [µm]	y [µm]		
1	-7542.31	14 -1806.	440	
3	-7752.92 -8206.76	22 -1970. 54 -1970.	.115 .115	
			Lo	ad Save

 Convex hull类似Bounding grid,通过 移动载物台到想拼接的图像边缘,单击Add ,添加边缘范围点。软件自动定义拼图范 围。不同于Bounding grid软件定义拼图范 围会更接近选择的位置点。
 "Start Experiment"。



1	Startup and Shutdown of the System
2	Acquiring Multi-Channel images
3	Z-stack image
4	Time Series image
5	Tile Scan
6	Airyscan Imaging
7	
8	









8、选择适合的滤片,如BP420-480代表蓝色波段,BP495-550代表绿色波段, BP 570-620 代表红色波段





9、次级分色镜(SBS)选择"Plate"





10、如果是多色荧光成像,可以选择多激光反射主分色镜(MBS), 如图中MBS可以同时反射488和561,减少MBS转换



11、不要更改针孔(Pinhole)大小, Airyscan Mode确认是"SR"模式

👻 🕿 Acquisition	Mode	🗸 Show all 🛛 📝
Objective	Plan-Apochromat 63x/1	I.4 Oil DIC M27
Scan Mode	Frame 🔻	
Frame Size	X 2000 🗘 🗙	* Y Y 2000 🗘
Line Step	1	Optimal
Speed	0	6 🗘 Max
Pixel Dwell	1.06 µsec Scan Time :	14.90 sec
Averaging		
Number	1 Bit	t Depth 8 Bit 🔹
Mode	Line 🔻 Di	rection <>
Method	Mean 🔻 Corr X	
	Corr Y	0.00
		Auto
HDR		
Airyscan Multitra	ck	ILEX Setup
😔 Scan Area		
(Image Size: 74	l.9 μm x 74.9 μm
	Pixel Size: 0.	04 μm
		0.0 ‡ C
		0.0 ¢ C
	• • •	0.0 \$ 0
	Zoom 🚽	
		Reset All

ZEISS

- 13、Frame Size选择"optimal"
- 14、Speed选择"Max"

12、Zoom不能小于1.8





15、随后可以进行二维图像拍摄(Snap)或者多维图像拍摄



III) ZEN 2			
File View Acquisition Maintain Macro Tools Window I	-Ip		
, , , , , , , , , , , , , , , , , , ,			
Locate Acquisition FCS Processing Maintain	Image 31.cl O Image 1 O	ZEN ZEN 2	
Apply	20	File View Acquisition Maintain Macro Tools Window He	
Method	Anystan		
Adjust Airyscan Processing Average	Spile		
Channel Alignment Color-coded projection		Locate Acquisition FCS Processing Maintain	
Copy Correlation			
Filter HDR- imaging	250	A Apply	
ICS Image calculator Incomposition	3		
	A Hiso	Method	
Input image I Select	Coloc	Adjust	
	ulua Profile	Airyscan Processing	
	FRET	Average	
Strength	Preview	Channel Alignment	
Run Batch	i Info	Color-coded projection	
Fleview		Correlation	
Input Preview			

17、按照"Processing → Airyscan Processing → Select → Apply" 进行图像处理。如果是z-stack图像,请选择"3D"

		Method Parameters	
		Input Image Image 1 3 Select	
		Output Image I_Airyscan Processing	
	Display 1 7 7 Z-Position 1 13 7 7 Zoom 0 1 125 1	Strength	-
		Run Batch	
	Channels Chick Carl Concest Single Channel Range Indicator Reuse S Crop Positio Stage		
PU 1% FreeHD 1.8TB eeRam 42GB			







18、处理后的图像,在Display处 选择Min/Max或者手动调节显示效 果

